



PREMIE HAL® S108.100

24-Week Preterm Newborn Skills Trainer

- Anatomically accurate oral cavity and airway
- Visible chest rise with assisted ventilation
- · Realistic chest compression and recoil
- Fontanelle, umbilical, brachial, and femoral pulses
- IV cannulation, umbilical artery, and catheterization

Critical Thinking and Skills Training Through Enriched Experiences

The Premie HAL® S108.100 is a lifelike, 24-week preterm patient simulator designed to facilitate the training of healthcare professionals in the areas of airway management, respiratory support, and resuscitation.











Suctioning

Bag ventilation

Active Resuscitation


Positive pressure ventilation

Realistic, Pre-Term Stabilization and Resuscitation Exercises

Premie HAL® S108.100 has a lifelike, anatomically accurate oral cavity and airway which allows participants to rehearse stabilization and resuscitation for preterm infants.

- Airway landmarks include: epiglottis, glottis, and vocal cords
- · Realistic chest recoil during CPR
- · Manual Pulses: Brachial, Femoral, Umbilical, and Fontanelle
- · Vascular access via hand, scalp, and umbilicus
- Supports:
 - » UVC/UAC sampling and infusion
 - » PICC line placement





FEATURES

- · Age: 24 Week preterm neonate
- Length 12.5" (31.75cm)
- Weight 1.3lbs (.59 kg)

AIRWAY

- · Lifelike, anatomically accurate oral cavity and airway
- Lifelike gums and appropriate sized tongue
- · Endotracheal intubation
- · Sellick maneuver
- · Nose and oral cavity suction

BREATHING

- · True-to-life lung compliance
- · Visible chest rise following recommended flow, PIP, and PEEP values
- Supports standard positive pressure ventilation devices including BVM, CPAP, and mechanical ventilators

CARDIAC

- · Pulses (manual)
 - » Brachial
 - » Femoral
 - » Umbilical
 - » Fontanelle
- · Realistic chest recoil during **CPR**

VASCULAR ACCESS

- IV cannulation
 - » Hand
 - » Scalp
- » UVC/UAC infusion and sampling
- PICC line placement
- Navel insert

GASTROINTESTINAL

- · Gastric distension
- · Patent esophagus
- NG/OG intubation
- · Gastric suctioning and feeding

True-to-Life Respiratory **Support Training**

Premie HAL® S108.100 has realistic lungs that respond to respiratory support like a real premature infant. This allows care providers to practice advanced, neonatal respiratory care in situ and interpret actual PV waveforms.

- True-to-life lung compliance
- · Visible chest rise following recommended flow. PIP. and PEEP values
- Supports standard positivepressure ventilation devices
 - » BVM
 - » CPAP
- » Mechanical ventilators

PREMIE HAL®

S108.100

Available skin tones





PACKAGE CONTENTS:

- · 24 week preterm neonate
- IV bag, stand, and tubing
- Filling kit
- · Drainage kit
- PICC line inserts
- · Squeeze bulb with tubing kit
- Navel insert
- Umbilical cords
- Replacement left IV arm
- · Baby blanket
- Diaper
- Mineral oil lubricant
- · Concentrated blood
- Medical tape
- · Canvas carrying bag

OPTION

GAUMARD VIRTUAL INSTRUMENTS* - PATIENT MONITOR SKILLS TRAINERS

REQUEST A QUOTE

14700 SW 136th Street Miami, Florida 33196-5691

P.O. Box 140098 Coral Gables, Florida 33114-0098

Toll Free USA Call 8:00 a.m. - 7:30 p.m. ET Monday - Friday 800.882.6655

Worldwide 305.971.3790

Fax 305.252.0755 24 hours

Website www.gaumard.com

Sales / Customer Service sales@gaumard.com

ABOUT GAUMARD®

Gaumard® is family owned and operated, and is the direct source for your health care education needs.

PLACE YOUR ORDER

By telephone, fax, or on-line at www.gaumard.com

WE ACCEPT

VISA, MasterCard, American Express

WARRANTY

Gaumard products are covered by a one-year limited warranty. Terms and conditions apply. Please visit www.gaumard.com for details.

TECHNICAL SUPPORT

Available 8:00 a.m. to 7:30 p.m. EST weekdays, or email us at support@gaumard.com Phone: 786.478.3838

ONLINE SUPPORT 24/7

Download articles, instructions, teaching tips. View instructional videos on demand.

REPAIR

Gaumard® offers repair services and parts. For more information visit www.gaumard.com

